

Product Warranty & Important Information

Your DITH™ EZ Sling is warranted against defects found in material and/or workmanship for six months when used in accordance with generally accepted medical use and in accordance with the instructions provided by DITH. Alteration to, misuse of, or failure to maintain use pursuant with the provided directions from DITH and your medical treating professional shall void this warranty. Normal wear and tear during the usage of a product is not considered a defect.

If it is determined that a return or replacement of the product is appropriate, instructions will be provided as to how to do so. If the product is unused and a request is made for return, it must be in the original packaging and in new resalable condition or no refund can be given. No merchandise is returnable without first obtaining company authorization.

What Patients are Saying:

"Helps decrease pain and pressure off the very painful area. Have seen wonderful results." ~ J.W.

"It helps ease pain. Numbness decreases with daily use. I wear it while watching t.v." ~ S.N.

"It helps me alot. It love it. It takes pressure off of my arm." ~ S.S.

"During times of prolonged pain or flare-ups, taking weight off shoulder aided in the relief of pain". ~ K.K.

Important Information

As with any medical device if you develop any new pain symptoms, you should contact your treating medical professional immediately. Symptoms such as pain, numbness, tingling, burning, chafing or irritation of the skin (sometimes found in the armpit due to lack of movement) and progressive stiffness are not normal. Doctor in the House™ makes every effort to ensure the accuracy and completeness of information provided with this product. Doctor in the House™ disclaims any liability or responsibility for injury or damage to persons or property which is incurred as a consequence, directly or indirectly, of the use and application of The DITH™ EZ Sling.

It is the users responsibility to know and follow local care protocol as provided by the medical advisors directing the system to which he or she belongs. Also, it is the reader's responsibility to stay informed of changes in the treatment of upper extremity injuries. Medicine is an ever changing science. In view of the possibility of human error or changes in medical science, neither the founders, employees or any other party who has been involved in the preparation or production of The DITH™ EZ Sling product, warrant that the design and information contained herein is in every respect accurate or complete and they are not responsible for any errors, omissions, defects or injuries from the results obtained from the use of information or The DITH™ EZ Sling.

Users are encouraged to confirm the information contained herein with other sources.

“Patent D514,224” “Patent D561,902”

1515 DeKalb Pike Suite 100A Blue Bell, PA 19422 • Phone 610-277-1990 x 151
visit us at www.docinthehouse.com



1515 Dekalb Pike
Suite 100A
Blue Bell, PA 19422
610-277-1990 x151

THE ULTIMATE ARM SLING, SHOULDER IMMOBILIZER

DITH™ EZ SLING

- SUPPORT FROM BELOW
- NO NECK STRAPS - eliminates neck pain
- ITS SMALL AND EASY TO USE



Recovering from shoulder injury or surgery has never been more comfortable and safe. Provides maximum support and sequential freedom in treatment for patients with Neck, Shoulder, Arm Injury.

THE WORLD'S SIMPLEST & MOST COMFORTABLE ARM SLING & SHOULDER IMMOBILIZER



- State-of-the-Art Device
- Greater Comfort
- More Convenient
- No Neck Straps
- Easy to use
- For All Body Types
- Used on right or left side

For questions, call our trained therapists who are here to help.

610-277-1990 x151 or 844-260-5162 (toll free)
or see our videos at www.docinthehouse.com

How to Apply your DITH™ EZ Sling in 4 Easy Steps

- 1** Remove the fanny pack from the plastic bag and clip the belt around your waist. Snug it up really tight for the best support. In this example, for right side use, rotate the bag around to the right front of your body (for left simply use left labeled straps).
- 2** As shown in the example to the right, place your right hand through the Right Biceps loop on the belt.
- 3** Place the arm board on top of the bag and slide your hand under the forearm strap (on top of the bag) from the outside, follow the arrow, toward your body and rest your forearm on top of the bag & arm board so that it is comfortable - fingers dangling over the edge.
- 4** Place the Right Wrist loop over the hand and arm board. Adjust the angle to your comfort and snug up the belt one last time really tight!



The optional abduction pillow attachment (housed in the bag) may be used to allow the shoulder to assume 40 degrees of abduction. This allows gentle lengthening of the capsular and brachial plexus nerve structures and just used 30 minutes once or twice a day helps prevent and or treat frozen shoulder.

You're all set!

