

Product Warranty & Important Information

Your Miracle Ankle from Doctor in the House™ is warranted against defects found in material and/or workmanship for six months when used in accordance with generally accepted medical use and in accordance with the instructions provided by Doctor in the House™. Alteration to, misuse of, or failure to maintain use pursuant with the provided directions from Doctor in the House™ and your medical treating professional shall void this warranty. Normal wear and tear during the usage of a product is not considered a defect.

If it is determined that a return or replacement of the product is appropriate, instructions will be provided as to how to do so. If the product is unused and a request is made for return, it must be in the original packaging and in new resalable condition or no refund can be given. No merchandise is returnable without first obtaining company authorization.

Important Information to Read

As with any medical device if you develop any new pain symptoms, you should contact your treating medical professional immediately. Symptoms such as pain, numbness, tingling, burning, chafing or irritation of the skin and progressive stiffness are not normal.

Doctor in the House™ makes every effort to ensure the accuracy and completeness of information provided with this product. Doctor in the House™ disclaims any liability or responsibility for injury or damage to persons or property which is incurred as a consequence, directly or indirectly, of the use and application of The Miracle Ankle from Doctor in the House™.

It is the users responsibility to know and follow local care protocol as provided by the medical advisors directing the system to which he or she belongs. Also, it is the reader's responsibility to stay informed of changes in the treatment of upper extremity injuries. Medicine is an ever changing science. In view of the possibility of human error or changes in medical science, neither the founders, employees or any other party who has been involved in the preparation or production of The Miracle Ankle from Doctor in the House™ product, warrant that the design and information contained herein is in every respect accurate or complete and they are not responsible for any errors, omissions, defects or injuries from the results obtained from the use of information for The Miracle Ankle from Doctor in the House™.

Users are encouraged to confirm the information contained herein with other sources.



1515 Dekalb Pike
Suite 100A
Blue Bell, PA 19422
610-277-1990 x 151

The Miracle Ankle from Doctor in the House™ *The Ultimate Brace for foot and ankle injury*

What would you say about a device that:

- Can treat fractures and ligament injury of the lower leg, ankle and foot
- Calms pain, eases stiffness from arthritis, tendinitis and helps fractures to heal
- Treats foot drop and ankle instability
- Was invented by an Orthopaedic Surgeon to help his patients heal without surgery
- That is clinically tested and proven to work

You'd call it a Miracle!



www.docinthehouse.com

What is the Miracle Ankle

- Ergonomically and Physiologically Correct - to actually **treat** foot drop, foot and ankle injuries
- Proper and adjustable position of the lower leg, foot and ankle
- Adjustable immobilization and fit

The Miracle Ankle from Doctor in the House is a unique Ankle Foot Orthosis developed by an orthopedic surgeon to help patients heal foot, ankle and lower leg injury without surgery and to help rehabilitate these patients, as well as post - operative patients to optimal health and lifestyle.

Indications For Use

Designed to treat talocalcaneal intraarticular injury, avascular necrosis of the talus, tibiotalar and talocalcaneal arthrosis and foot drop.

The device can be adjusted to limit motion in the case of fracture healing and/or articular defects to enhance healing in optimal position and avoid progressive trauma to intraarticular fractures.



Foot Drop

- Chronic medial or lateral ankle
- Instability Arthritis and DJD of the ankle joint Fractures and Dislocations
- Strain and Sprain of the ankle
- Lower leg and ankle Muscle Injury and Tendinitis

Application of the Miracle Ankle

How to apply lower leg portion as healing progresses

Step 1: Unfasten top, middle and lower cuffs.

Step 2: Place the foot and ankle in the brace, as you would a boot. Position the heel so it is against the back of the soft liner.

Step 3: Close the soft liner and fasten the Velcro straps snugly to comfort.




Make sure the brace still fits snug and is centered over each side of the ankle.

Posterior Panel Application

For added support, the posterior panel can be attached to the ankle brace.

Step 1: The large posterior panel affords maximum support for ankle injury and foot drop. Secure the panel via the multiple Velcro attachment points to fit your leg length and replace the inner liner.



Step 2: Secure the brace as above (picture on the right) and additionally secure the upper straps.

YOU'RE ALL SET!

***Notice:** This device is intended for use only under the direction of a medical professional. If you experience any pain, swelling, sensation changes, or any unusual reactions while using the product, consult your medical professional immediately.

Cleaning instruction: Remove soft liner from brace and hand wash in cold water with mild detergent. Wipe the remainder of the brace with cold water. Allow all parts to air dry before re-assembly. Do not dry clean, iron, or bleach any part of this brace, as it could cause damage and void the manufacturer's warranty.