

## Product Warranty & Important Information

Your DITH™ Carpal Cure is warranted against defects found in material and/or workmanship for six months when used in accordance with generally accepted medical use and in accordance with the instructions provided by DITH. Alteration to, misuse of, or failure to maintain use pursuant with the provided directions from DITH and your medical treating professional shall void this warranty. Normal wear and tear during the usage of a product is not considered a defect.

If it is determined that a return or replacement of the product is appropriate, instructions will be provided as to how to do so. If the product is unused and a request is made for return, it must be in the original packaging and in new resalable condition or no refund can be given. No merchandise is returnable without first obtaining company authorization.

### Important Information to Read

As with any medical device if you develop any new pain symptoms, you should contact your treating medical professional immediately. Symptoms such as pain, numbness, tingling, burning, chafing or irritation of the skin and progressive stiffness are not normal.

Doctor in the House™ makes every effort to ensure the accuracy and completeness of information provided with this product. Doctor in the House™ disclaims any liability or responsibility for injury or damage to persons or property which is incurred as a consequence, directly or indirectly, of the use and application of The DITH™ Carpal Cure.

It is the users responsibility to know and follow local care protocol as provided by the medical advisors directing the system to which he or she belongs. Also, it is the reader's responsibility to stay informed of changes in the treatment of upper extremity injuries. Medicine is an ever changing science. In view of the possibility of human error or changes in medical science, neither the founders, employees or any other party who has been involved in the preparation or production of The DITH™ Carpal Cure product warrant that the design and information contained herein is in every respect accurate or complete and they are not responsible for any errors, omissions, defects or injuries from the results obtained from the use of information or The DITH™ Carpal Cure.

Users are encouraged to confirm the information contained herein with other sources.

Patents 9,173,764. 9,895,251. 0779,075. Patents Pending



1515 Dekalb Pike  
Suite 100A  
Blue Bell, PA 19422  
610-277-1990 x 151

### The DITH™ Carpal Cure

*The Ultimate Brace for Hand & Wrist Injury*

#### What would you say about a device that:

- Can treat fractures of the hand and wrist, tendinitis and carpal tunnel
- Calms pain, eases stiffness from arthritis, tendinitis and helps fractures to heal
- Helps patients to heal happily with lessened or eliminated discomfort
- Was invented by an Orthopaedic Hand Surgeon to help his patients heal without surgery
- That is clinically tested and proven to work

#### You'd call it a Carpal Tunnel Cure!



[www.docinthehouse.com](http://www.docinthehouse.com)

## What is the DITH™ Carpal Cure

- Customizable - microwave 1 minute and custom molds to fit **you** like a glove
- Ergonomically and Physiologically Correct - to actually **treat** hand, wrist and forearm injuries
- Proper neutral position of the wrist
- Freedom of motion for activities

The purposeful addition of multiple treatment and support layers allows trapping of the body's natural thermal energy and muscular contractions. It provides soothing heat, massage and comfort to the forearm, hand and wrist structures while providing protection and stability to the injured tissues through the heatable, moldable, removable and changeable internal thermoplastic struts.

### Indications For Use

- Carpal Tunnel Syndrome
- Repetitive Strain Injury
- Fractures about the wrist, hand, fingers and forearm
- Arthritis
- Tendinitis or Tenosynovitis about the hand, wrist and forearm
- Strain and Sprain about the wrist and hand
- Lateral Epicondylitis (Tennis Elbow), Medial Epicondylitis
- Extensor and Flexor Tendinitis about the hand & wrist


### How to use your Carpal Cure

Congratulations on your first step toward self healing. The brace feels good and helps you to heal. Follow a few simple steps and you are on your way.




## How to use your Carpal Cure

### Using your Carpal Cure is as easy as 1, 2, 3



#### Step 1:

Remove the Carpal Cure brace from packaging and use velcro straps to adjust to your injured area.



#### Step 2:


Place the Carpal Cure brace in microwave oven for 60 seconds. Make sure temperature is not too hot before slipping your thumb and hand in brace. Tighten the straps.



#### Step 3:

Rest your hand and forearm palm down on a table for 60 seconds to custom mold to your hand and wrist.

**YOU'RE ALL SET!**



### Optional Night Rest

Simply slide out piece A, insert the long end of piece D, secure the flap and slide in your hand and fingers to rest on the night rest!



Shown with optional night rest for arthritis

Please refer to included DVD and your medical professional for hinge and dynamic traction.