

Product Warranty & Important Information

Read and Follow the Molding and Heating instructions that are included

Your DITH™ Miracle Splint Dynamic is warranted against defects found in material and/or workmanship for six months when used in accordance with generally accepted medical use and in accordance with the instructions provided by DITH. Alteration to, misuse of, or failure to maintain use pursuant with the provided directions from DITH™ and your medical treating professional shall void this warranty. Normal wear and tear during the usage of a product is not considered a defect.

If it is determined that a return or replacement of the product is appropriate, instructions will be provided as to how to do so. If the product is unused and a request is made for return, it must be in the original packaging and in new resalable condition or no refund can be given. No merchandise is returnable without first obtaining company authorization.

Important Information to Read

As with any medical device, if you develop any new pain symptoms, you should contact your treating medical professional immediately. Symptoms such as pain, numbness, tingling, burning, chafing or irritation of the skin and progressive stiffness are not normal.

Doctor in the House™ makes every effort to ensure the accuracy and completeness of information provided with this product. Doctor in the House™ disclaims any liability or responsibility for injury or damage to persons or property which is incurred as a consequence, directly or indirectly, of the use and application of The DITH™ Miracle Splint Dynamic.

It is the users responsibility to know and follow local care protocol as provided by the medical advisors directing the system to which he or she belongs. Also, it is the reader's responsibility to stay informed of changes in the treatment of upper extremity injuries. Medicine is an ever changing science. In view of the possibility of human error or changes in medical science, neither the founders, employees or any other party who has been involved in the preparation or production of The DITH™ Miracle Splint Dynamic product warrant that the design and information contained herein is in every respect accurate or complete and they are not responsible for any errors, omissions, defects or injuries from the results obtained from the use of information or The DITH™ Miracle Splint Dynamic.

Users are encouraged to confirm the information contained herein with other sources.

Patents 9,173,764 9,895,251 D779,075 D824527S



1515 Dekalb Pike
Suite 100A
Blue Bell, PA 19422
610-277-1990 x 151

THE MIRACLE Splint Dynamic

The Answer to Every Hand/Wrist Problem:

- Helps Heal Carpal Tunnel, Tendon Injury, Tenosynovitis, Wrist Injury, Arthritis and Tendinitis
- Dynamic Splinting for Zone 2, 3, 4 Tendon Injury Fractures, Tendinitis, Bone and Ligament Injury
- Elastic Motion and Hinge Stabilization of Wrist Fractures and Digital Arthritis
- Invented and tested by an Orthopaedic Hand Surgeon to help his patients heal without surgery
- Locking hinge to tx Lunatotriquetral Ligament

The Ultimate Treatment Miracle!



docinthehouse.com HCPCS L3915 & L3916

What is the DITH™ Miracle Splint Dynamic?

- Dynamic Traction, Locking Hinge and Customizable
- Treats injury with heat, massage & protected motion - actually **treats** hand, wrist, and forearm injuries
- Wrist Motion Stabilization as needed with locking hinge
- Dynamic Healing and Stabilization

Mold and Customize

Your doctor, therapist or orthotist will customize and mold your device as well as the struts, hinge and dynamic traction to optimally treat your specific injury or medical problem.

Institute Dynamic Traction. Hinge institutes varied locked and progressive motion. Customized to fit with molding, trimming and molding of the traction tower and locking hinge to address the specific injury and medical indications.

Indications For Use:

- Lunatotriquetral Ligament Injury and Fractures
- Extensor and Flexor Tendon Lacerations & Post Op
- Fractures about the wrist, hand, fingers and forearm
- Tendinitis, Arthritis or Tenosynovitis about the hand, wrist, and forearm
- Carpal Tunnel, Trigger Digits and RSI
- Lateral Epicondylitis (tennis elbow), Medial Epicondylitis
- Zone 2,3,4 Flexor and Extensor Tendon Injury and Post op

How to use your Miracle Splint Dynamic?

Using your Miracle Splint is as easy as 1, 2, 3



Step 1:

Remove the brace from packaging and use velcro straps to adjust to your injured area.



Step 2:

Place the Miracle Splint in microwave oven for 60 seconds. Make sure temperature is not too hot before slipping your thumb and hand in brace. Tighten the straps.


THE BACK SHOULD BE OPEN



Step 3:

Rest your hand and forearm palm down on a table for 60 seconds to custom mold to your hand and wrist.

YOU'RE ALL SET!



SUPPORT IS FROM THE PALM SIDE



Optional Night Rest

Simply slide out piece A, insert the long end of piece D, secure the flap and slide in your hand and fingers to rest on the night rest!



Shown with optional night rest for arthritis.

Please refer to customizing instructions and your medical professional for hinge and dynamic traction.